


Omni-Flex™

Urology Retractor System



Your Patient.
Your Approach.

Flexibility

Interchangeable, Removable Arms

Stable retraction that lets you perform procedures the way you prefer. Strong enough for your largest patients. Nimble enough for small-site work. Always easy to use.

Use one arm. Use two arms. Take one arm away or use one of each. Choose the frame arm—before or during surgery—that helps maximize your surgical field.

Ease-of-Use


Maximize Exposure—Minimize Obstructions

Exchange an arm. Remove an arm. Open the field. The Omni-Flex System easily exposes difficult anatomy no matter what your procedure demands.

Strength

Strength for Stability with Every Patient

Any patient size. Any incision size. Always stable. The Omni-Flex System is designed to remain stable, so you always get the exposure you need.



Retract from every direction by simply adding our hinged extension arm.



flexibility

ease-of-use



Hands-free exposure during perineal and



SideSnap™ Clamps for Faster Setup & Mid-Surgery Changes

Snap the clamp on the frame. Snap the blade in place. Our SideSnap™ Clamps allow for quicker, simpler setup. With added strength, they enhance the flexibility of the Omni-Flex System.



Quick Blade Attachment

Attach in one sweeping motion with a simple snap from the side.



Stack Blades to Maximize Incision


Stack a longer blade over a shorter blade for increased tissue holding power and increased visibility.

of-use
strength



minimally invasive abdominal approaches.

Product Specifications



Omni-Flex™ Urology Retractor Systems

U0550 Expanded Urology Retractor System
 U0500 Urology Retractor System
 U0100 Pelvic Retractor System

	U0550	U0500	U0100	Description	Part Number
A	1	1	1	Sterile Field Post	10244
B	1	1	1	Support Arm	11050
C	2	2	--	Standard WISHBONE® Frame Arm	10269
D	2	2	2	Small WISHBONE® Frame Arm	10270
E	8	6	5	SideSnap™ Clamp	11030
F	1	--	--	Hinged Extension Arm, 19" (48.3cm)	11121
G	--	--	1	Extension Arm, 20" (50.8cm)	11070
H	--	--	1	Anterior Swivel Retractor, 1½" x 3" (3.8cm x 7.6cm)	4977
I	1	1	--	Deaver Swivel Retractor, 2¾" x 4½" (7cm x 11.4cm)	3784
J	1	1	--	Deaver Swivel Retractor, 1½" x 6" (3.8cm x 15.2cm)	3854
K	1	1	--	Fence Swivel Retractor, 4" x 6" (10.2cm x 15.2cm)	3772
L	1	1	--	Harrington Sweetheart Swivel Retractor, 2" x 5" (5cm x 12.7cm)	3794
M	1	1	--	Harrington Sweetheart Swivel Retractor, 2½" x 7" (6.4cm x 17.8cm)	3796
N	--	--	4	Malleable Swivel Retractor, ¾" x 2" (1.9cm x 5cm)	4964
O	--	--	4	Malleable Swivel Retractor, ¾" x 3" (1.9cm x 7.6cm)	4965
P	--	--	4	Malleable Swivel Retractor, ¾" x 4" (1.9cm x 10.2cm)	4966
Q	1	1	--	Malleable Slotted Swivel Retractor, 2" x 6½" (5cm x 16.5cm)	3788
R	1	1	--	Malleable Slotted Swivel Retractor, 3" x 6½" (7.6cm x 16.5cm)	3790
S	1	1	--	Marshall Bladder Malleable Swivel Retractor, 3" x 6" (7.6cm x 15.2cm)	3792
T	2	2	--	Mayo Swivel Retractor, 2" x 2" (5cm x 5cm)	3883
U	2	2	--	Mayo Swivel Retractor, 2¾" x 2" (7cm x 5cm)	3335
V	--	--	1	Posterior Swivel Retractor, 1¾" x 3" (4.5cm x 7.6cm)	4975
W	--	--	1	Posterior Swivel Retractor, 1¾" x 4" (4.5cm x 10.2cm)	4976
X	1	1	--	Splanchnic Swivel Retractor, 1½" x 4" (3.8cm x 10.2cm)	3798



4849 White Bear Parkway, St. Paul, MN 55110 USA
 800-367-8657 • 651-287-4300 • fax: 651-287-4400
 www.omni-tract.com

© 2006 Omni-Tract Surgical
 All Rights Reserved.

U.S. Patent Nos. 4,355,631; 4,617,916; 4,718,151; 4,930,932; 4,949,707;
 5,020,195; 5,242,240; 5,380,338; 5,400,772; 5,688,223; 5,727,899;
 5,792,046; 5,882,298; 5,899,627; 6,264,396; 6,315,718; 6,663,563,
 D343,235; Other Patents Pending. Omni-Flex and SideSnap are trademarks,
 Omni-Tract and WISHBONE are registered trademarks of Minnesota Scientific, Inc.



3-06/40-6155/1000