



GhosTract™

Radiolucent Retractor Blade Kit




Now you see them...



with GhosTract™, you don't!

- Perfect for anterior lumbar spine procedures
- Virtually artifact free for imaging
- Cost effective re-usable aluminum construction
- Compatible with both "ring" and "wishbone" style retractor systems

Product Specifications



5605 GhosTract™ kit, square blade holder for ring retractors
5610 GhosTract™ kit, round blade holder for wishbone retractors
Includes:

Qty			
A.	1	Deaver, radiolucent 2 3/4" x 4 1/2" (7cm x 11.4cm)	11300
B.	1	Splanchnic, radiolucent 2" x 4 1/2" (5cm x 11.4cm)	11308
C.	1	Splanchnic, radiolucent 2" x 6" (5cm x 15.2cm)	11309
D.	2	Malleable, radiolucent 2" x 6 1/2" (5cm x 16.5cm)	11301
E.	2	Renal, radiolucent 1" x 8" (2.5cm x 20.3cm)	11306
F.	2	Renal, radiolucent 1" x 6" (2.5cm x 15.2cm)	11305
G.	2	Renal, radiolucent 1" x 4" (2.5cm x 10.2cm)	11304
H.	2	Mayo/Balfour, radiolucent 2 3/4" x 2" (7cm x 5cm)	11303
I.	2	Mayo, radiolucent 3 1/2" x 2" (9cm x 5cm)	11302
J.	6	Blade Holder, square	11340
		or	
K.	6	Blade Holder, round	11360
L.	1	Sterilization tray	11355



4849 White Bear Parkway, St. Paul, MN 55110 USA
 800-367-8657 • 651-287-4300 • fax: 651-287-4400
 www.omni-tract.com

© 2007 Omni-Tract Surgical
 All Rights Reserved.

U.S. Patent Nos. 4,355,631; 4,617,916; 4,718,151; 4,930,932; 4,949,707;
 5,020,195; 5,242,240; 5,380,338; 5,400,772; 5,688,223; 5,727,899; 5,792,046;
 5,882,298; 5,899,627; 6,264,396; 6,315,718; 6,663,563, D343,235; Other
 Patents Pending. Omni-Flex and SideSnap are trademarks; Omni-Tract and
 WISHBONE are registered trademarks of Minnesota Scientific, Inc.



6/07/40-9005/2500